

4D-Nucleofector® 96-well Unit Guideline for Bacteria Transformation

1. Kit for Bacteria Transformation:

We do not offer a kit specialized for bacteria transformation. Please contact your local sales representative or our Scientific Support Team for a special Sales Agreement to access 96-well Nucleocuvette® Plates as custom product.

2. Protocol

Programs:

- Start your respective Nucleofector® System and create or upload the experimental parameter file (see system manual for details)
- This special program definition file contains the standard 96-well programs for eukaryotic cells plus **17 different programs for bacteria** which can be accessed most easily via the “cell type code” drop down:

Table 1:

Cell type code	Program Code
X_Bacteria 1	GA-100
X_Bacteria 2	GB-100
X_Bacteria 3	GC-100
X_Bacteria 4	GD-100
X_Bacteria 5	GE-100
X_Bacteria 6	GF-100
X_Bacteria 7	GG-100
X_Bacteria 8	GH-100
X_Bacteria 9	GI-100
X_Bacteria 10	GJ-100
X_Bacteria 11	GK-100
X_Bacteria 12	GL-100
X_Bacteria 13	GM-100
X_Bacteria 14	GN-100
X_Bacteria 15	GO-100
X_Bacteria 16	GP-100
X_Bacteria 17	GQ-100

We recommend following optimization strategies:

- If you work with an *E.coli* strain please try programs GM-100, GN-100 and GO-100
- For any other strain we recommend an initial screen using all 17 programs
- Programs GO-100, GP-100 and GQ-100 may lead to error 8, depending on the sample background

Bacteria stocks:

- The Nucleofection® Procedure (see below) should be the same for any strain. Nevertheless the competence (e.g. # of colonies/μg pUC19) for electroporation as well as for chemical/heat shock transformation differs significantly for different strains
- Bacteria stocks should not contain more than 4–5% glycerol as pipetting of small volumes with higher glycerol amounts could lead to use of imprecise volumes.

DNA stock:

Ideally, the plasmid DNA should be dissolved in water. Use of any buffers may lead to arcing or splashing.

Nucleofection® Procedure:

1. Thaw electro-competent bacteria on ice
2. Pre-warm SOC medium to 37°C
3. Optionally: Pre-cool 96-well Nucleocuvette® Plate on ice for at least 15 minutes. With pre-cooled vessels you may increase transformation efficiencies
4. Dispense desired amount of plasmid DNA (e.g. 10 – 50 pg pUC19; in 2 – 5 μl of H₂O) in each well of the 96-well Nucleocuvette® Plate
5. Add 20 μL of bacteria suspension to each well of the 96-well Nucleocuvette® Plate
6. Carefully remove any ice-water from the plate

7. Place the 96-well Nucleocuvette® Plate into the 4D-Nucleofector® 96-well Unit and start the Nucleofection® Process as described in the system manual
8. Post Nucleofection®, remove the 96-well Nucleocuvette® Plate and add 200 µL SOC medium (or your established medium) to each well
9. Transfer 200 µL bacteria suspension to a 10 ml bacteria tube containing 300 µL SOC (total volume 500 µL)
10. Incubate at 37°C and 225 rpm for 1 h
11. Plate 25 µL or 50 µL of the bacteria suspension of each tube on an agar plate (containing appropriate antibiotics and 40 µg/ml X-gal and 100 µM IPTG for blue/white screen (if desired))
12. Incubate at 37°C for additional 24 h

Additional Information

For additional information and an up-to-date list of Nucleofector® References, please visit the Lonza Knowledge Center:

<https://knowledge.lonza.com>

For more technical assistance, contact our Scientific Support Team:

USA /Canada

Phone: 800 521 0390 (toll-free)

Fax: 301 845 8338

Email: scientific.support@lonza.com

Europe and Rest of World

Phone: +49 221 99199 400

Fax: +49 221 99199 499

Email: scientific.support.eu@lonza.com

Lonza Cologne GmbH – 50829 Cologne Germany

4D-Nucleofector® is a (registered) trademark of Lonza

For research use only. Not for use in diagnostic procedures. The 4D-Nucleofector® is covered by patent and/or patent pending rights owned by the Lonza Group Ltd or its affiliates. All trademarks belong to Lonza or its affiliates, and are registered in the USA, EU and/or CH, or belong to their respective third-party owners and are only being used for informational purposes. All third-party copyrights have been reproduced with permission from their owners. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. Any user must make his own determination and satisfy himself that the products supplied by Lonza Group Ltd or its affiliates and the information and recommendations given by Lonza Group Ltd or its affiliates are (i) suitable for intended process or purpose, (ii) in compliance with all applicable laws, including all environmental, health and safety regulations, and (iii) will not infringe any third party's intellectual property rights. The user bears the sole responsibility for determining the existence of any such third-party rights, as well as obtaining any necessary licenses and approvals. For more details: <http://www.lonza.com/legal>.

©2024 Lonza. All rights reserved. CD-OP034 12/24