

MACHEREY-NAGEL

# Universal reagent for RNA isolation

Bioanalysis




















NucleoZOL – Do it right!

- No chloroform, no phase separation, easy handling
- High RNA yield and purity from any sample material
- Combination with proven NucleoSpin® technology available

# NucleoZOL – The universal reagent for RNA isolation

## Procedure

NucleoZOL		Competitor Zol	
	Sample homogenization		Sample homogenization
	Addition of water + Non-toxic		Addition of chloroform - Toxic
	Centrifugation at RT Precipitation of contaminants		Centrifugation at 4 °C Liquid-liquid phase separation
	Aspiration of whole supernatant with RNA + Easy sampling of RNA + DNA and proteins remain in pellet		Aspiration of aqueous phase with RNA - Inconvenient pipetting for phase separation - Risk of carry-over of interphase / polar phase - Refrigerated centrifuge necessary
Standard procedure		Procedure with NucleoSpin® RNA Set for NucleoZOL	
	Centrifugation at RT Precipitation of RNA		Binding of RNA + Short centrifugation steps (2 x 30 sec)
	Washing of RNA RNA in pellet		Washing of RNA
	Resuspension of RNA + No drying of RNA necessary + Quick and easy		Elution of RNA + Proven NucleoSpin® purity + Standardized procedure
			Centrifugation at 4 °C Precipitation of RNA
			Washing of RNA RNA in pellet
			Drying and resuspension of RNA - Time-consuming drying of RNA

### Summary of advantages of NucleoZOL compared to competitor Zol

The table indicates the advantages of NucleoZOL procedure compared to typical competitor Zol products. The NucleoZOL procedure offers an easy and much more convenient liquid handling. A laborious chloroform two-phase separation is not necessary. Additionally, NucleoZOL can be combined with the proven NucleoSpin® technology.


	NucleoZOL	Common competitor products
Procedure	One phase procedure minimizes the risk of contamination by carry-over	Phase separation leads to risk of DNA / protein / phenol carry-over or sample loss
Removal of DNA and proteins	Precipitation of DNA and proteins results in minor risk of contamination	Contamination possible due to difficult phase separation
Solubilization of RNA	No drying required	Time-consuming drying step required
miRNA isolation	Protocol for fractionation of small and large RNA	No protocol for selective miRNA isolation available
Handling	All steps are performed at room temperature	Necessity of refrigerated centrifuge

# NucleoZOL – The universal reagent for RNA isolation

## Product at a glance


Specifications																					
Technology	One-phase separation																				
Sample material/mL reagent (scalable)	< 1 x 10 <sup>7</sup> cultured cells, bacteria, and yeast, < 100 mg human / animal / plant tissue, < 0.4 mL (viral) fluids																				
Fragment size	Small RNA = 10–200 nt, large RNA > 200 nt																				
Typical yield	<table border="0"> <tr> <td>Total RNA:</td> <td></td> <td>Large RNA:</td> <td></td> </tr> <tr> <td>Liver:</td> <td>6–8 µg/mg tissue</td> <td>Liver:</td> <td>5–7 µg/mg tissue</td> </tr> <tr> <td>Kidney, spleen:</td> <td>3–4 µg/mg tissue</td> <td>Kidney, spleen:</td> <td>3–4 µg/mg tissue</td> </tr> <tr> <td>Muscle, brain, lung:</td> <td>0.5–1.5 µg/mg tissue</td> <td>Muscle, brain, lung:</td> <td>0.5–1.5 µg/mg tissue</td> </tr> <tr> <td>Cultured cells:</td> <td>4–10 µg/10<sup>6</sup> cells</td> <td>Cultured cells:</td> <td>3–8 µg/10<sup>6</sup> cells</td> </tr> </table>	Total RNA:		Large RNA:		Liver:	6–8 µg/mg tissue	Liver:	5–7 µg/mg tissue	Kidney, spleen:	3–4 µg/mg tissue	Kidney, spleen:	3–4 µg/mg tissue	Muscle, brain, lung:	0.5–1.5 µg/mg tissue	Muscle, brain, lung:	0.5–1.5 µg/mg tissue	Cultured cells:	4–10 µg/10 <sup>6</sup> cells	Cultured cells:	3–8 µg/10 <sup>6</sup> cells
Total RNA:		Large RNA:																			
Liver:	6–8 µg/mg tissue	Liver:	5–7 µg/mg tissue																		
Kidney, spleen:	3–4 µg/mg tissue	Kidney, spleen:	3–4 µg/mg tissue																		
Muscle, brain, lung:	0.5–1.5 µg/mg tissue	Muscle, brain, lung:	0.5–1.5 µg/mg tissue																		
Cultured cells:	4–10 µg/10 <sup>6</sup> cells	Cultured cells:	3–8 µg/10 <sup>6</sup> cells																		
A <sub>260</sub> /A <sub>280</sub>	1.8–2.1																				
Typical RIN	> 9																				
Elution volume	Flexible																				
Preparation time	< 1 h																				

## Application data




### High RNA quality from fibrous tissue

Total RNA was isolated with NucleoZOL from 60 mg mouse heart tissue. RNA was analyzed on an Agilent Bioanalyzer. The RIN of 9 indicates perfect RNA quality.




### Low DNA contamination with NucleoZOL

DNA contamination was quantified by qPCR (competitor Zol=100%). Compared to a standard two-phase extraction, only a minimum of DNA is carried over during purification with NucleoZOL.



### Market-leading RNA yields

RNA was extracted from different starting materials. RNA was quantified by qRT-PCR and relative yields were calculated (NucleoZOL=100%). RNA isolation with NucleoZOL results in similar or better RNA yields compared to standard two-phase extraction methods (competitor Zol).



### High superior RNA recovery and easy handling

Total RNA was isolated with NucleoZOL combined with mini spin columns (NucleoSpin® RNA Set for NucleoZOL) or with a competitor product, that combines a liquid extraction with a column-based method (competitor Z). Total RNA yields were quantified using an Agilent 2100 Bioanalyzer™.

Product	Preps/Pack of	REF
NucleoZOL	200 mL	740404.200
NucleoSpin® RNA Set for NucleoZOL	10/50	740406.10/.50

# NucleoZOL – The universal reagent for RNA isolation

## Do it right!

During the years, MN Bioanalysis has gained vast experience in nucleic acid purification and thus developed into highly skilled RNA experts. The MN research and the technical support team are current in all RNA applications. This has enabled MN to provide high-value RNA purification products in view of the expansive research applications.

We invite you to take advantage of the excellent RNA purification products as well as our team of scientific experts from MACHERY-NAGEL.



- Qualified
- Customer-focused
- Reliable

Our friendly team is looking forward to give you professional advices to our wide range of products!

+49 24 21 969-0  
tech-bio@mn-net.com

## Ordering information

Product	Preps / Pack of	REF
<b>RNA from cells and tissue</b>		
NucleoSpin® RNA Plus	10 / 50 / 250	740984.10 / .50 / .250
NucleoSpin® RNA	10 / 50 / 250	740955.10 / .50 / .250
NucleoZOL	200 mL	740404.200
NucleoSpin® RNA Set for NucleoZOL	10 / 50	740406.10 / .50
NucleoSpin® RNA XS	10 / 50 / 250	740902.10 / .50 / .250
NucleoSpin® RNA Midi	20	740962.20
NucleoSpin® 8 RNA	12 x 8 / 60 x 8	740698 / .5
NucleoSpin® 8 RNA Core Kit	48 x 8	740465.4
NucleoSpin® 96 RNA	2 x 96 / 4 x 96 / 24 x 96	740709.2 / .4 / .24
NucleoSpin® 96 RNA Core Kit	4 x 96	740466.4
NucleoMag® 96 RNA	1 x 96 / 4 x 96	744350.1 / .4
<b>MicroRNA</b>		
NucleoSpin® miRNA	10 / 50 / 250	740971.10 / .50 / .250
NucleoSpin® miRNA Plasma	10 / 50 / 250	740981.10 / .50 / .250
Exosome Precipitation Solution (Serum/Plasma)*	2 mL / 12 mL / 60 mL	740398.2 / .12 / .60
Exosome Precipitation Solution (Urine)*	12 mL / 50 mL / 250 mL	740399.12 / .50 / .250
<b>RNA, DNA, and protein</b>		
NucleoSpin® TriPrep	10 / 50 / 250	740966.10 / .50 / .250
NucleoSpin® RNA/Protein	10 / 50 / 250	740933.10 / .50 / .250
NucleoSpin® RNA/DNA Buffer Set	100	740944
<b>RNA from blood</b>		
NucleoSpin® RNA Blood	10 / 50 / 250	740200.10 / .50
NucleoSpin® RNA Blood Midi	20	740210.20
NucleoSpin® 8 RNA Blood	12 x 8 / 60 x 8	740220 / .5
NucleoSpin® 96 RNA Blood	2 x 96 / 4 x 96	740225.2 / .4
<b>Small and large RNA from FFPE samples</b>		
NucleoSpin® totalRNA FFPE	10 / 50 / 250	740982.10 / .50 / .250
NucleoSpin® totalRNA FFPE XS	10 / 50 / 250	740969.10 / .50 / .250
<b>RNA from plant</b>		
NucleoSpin® RNA Plant	10 / 50 / 250	740949.10 / .50 / .250
<b>Poly(A) mRNA isolation from total RNA</b>		
NucleoTrap® mRNA Mini	12	740655
NucleoTrap® mRNA Midi	12	740656

\*This product is not available in the USA

www.mn-net.com

# MACHERY-NAGEL



MACHERY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany

DE / International:

Tel.: +49 24 21 969-0

Fax: +49 24 21 969-199

E-mail: info@mn-net.com

CH:

Tel.: +41 62 388 55 00

Fax: +41 62 388 55 05

E-mail: sales-ch@mn-net.com

FR:

Tel.: +33 388 68 22 68

Fax: +33 388 51 76 88

E-mail: sales-fr@mn-net.com

US:

Tel.: +1 484 821 0984

Fax: +1 484 821 1272

E-mail: sales-us@mn-net.com

