

# scientifix

PRESS RELEASE.

## DYNAMICA pH MASTER, LOW COST HIGH PERFORMANCE pH METER

The New pHMaster range of pH meters from Dynamica and available from Scientifix is designed for the everyday lab use and also specialized applications in biological and life sciences.

The small footprint belies the advanced features of the pHMaster.

An easy to use keypad layout makes the systems easy to use and calibrate and the robust casing withstands the rigors of the laboratory and makes the unit ideal for teaching labs.

All models feature a user selectable 1, 2 or 3 points calibration with an automatic buffer recognition function and the three decimal place readout makes the unit flexible enough for even the most rigorous research uses.

The pHMaster BIO is designed with biological and life science applications in mind and is packaged with a combination electrode specific for TRIS buffer preparation thus overcoming the issues of inaccuracy obtained with general pH probes.

There is also an optional micro pH electrode with a 4.5mm diameter for small volume applications in 1.5/2ml micro tubes.

For general laboratory applications, the pHMaster LAB is the ideal choice as it is packaged with a convenient 3-in-1 pH/ATC electrode which continuously compensates for any temperature fluctuations. Furthermore a variety of different electrodes are available for many applications such as environmental water, food, solid or semi-solid samples.

The optional smart stirrer provides the homogeneous conditions for more accurate pH measurements.

All this at a very affordable price!

For more information concerning these innovative products please visit the Dynamica website

<http://www.dynamica-eu.com/cpzs5.asp?sortsid=837>

For more information please contact

Ken Baker at Scientifix on 1800 007 900 or 03 8540 5900

Email [ken@scientifix.com.au](mailto:ken@scientifix.com.au)

