

scientifix

PRESS RELEASE.

INCUBATORS IN A VARIETY OF SIZES AND FORMATS FOR ALL PURPOSES.

A wide range of Incubators designed for specific purposes are available from Scientifix Pty Ltd.

Scientifix as you may be aware is a specialist high tech equipment supplier and prides itself on the knowledge base of our skilled staff to give valuable advice on the selection of equipment. In science there is never one size that fits all and in incubation this is certainly true.

Scientifix source incubators from quality suppliers including Froilabo, N Biotek, Heal Force and S.E.M Pty Ltd.

GENERAL INCUBATORS

Froilabo produce a wide range of gravity convection, forced air and refrigerated incubators in sizes from 60 litres to 240 litres and provide an excellent, general purpose incubator at highly competitive prices. For larger sizes in general incubation we also supply the Australian made S.E.M units from Adelaide. These are ideal for routine microbiology, drying and teaching applications.

For more information visit

<http://www.froilabo.com/etuve-incubatoare-bacteriological-incubators.php>

and

www.semsa.com.au

CO2 AND / OR SHAKING INCUBATORS

For Shaking incubators and simple CO2 incubators, the N Biotek range is an exceptional product. These incubators are available in sizes from 60 litres up to 800litres with or without CO2 and shaking, and offer six sided heating for maximum temperature uniformity and rapid heat recovery. They have a small footprint and can be stacked and the high quality manufacturing ensures that they will give many years of trouble free service.

Whether your need is for one small incubator in the corner or several high volume multi layered stacked, stirring incubators for tissue culture, N Biotek will please you with their value for money.

For more information visit

http://www.n-biotek.com/?module=Board&action=SiteBoard&iBrdNo=1&sMode=SELECT_FORM

CO₂, WATERJACKETED AND TRIGAS INCUBATORS.

HealForce manufacture a range of CO₂, water jacketed and Tri gas incubators for the most demanding of incubation requirements. If you have a need for an incubator which includes 90degree moist heat disinfection the HF90 or HF 240 models are the ideal choice. Their easy clean design ensures that contamination control is simple and design features such as the direct heating and efficient humidity control mean that conditions within the incubator are controlled efficiently and a series of alarms alert the operator to any conditions which may compromise the incubator contents.

For more information visit the following web link

<http://www.medicaldeviceschina.com/2-4-incubator.html>

All incubation products are fully supported by Scientifix and we can offer a full contamination control service as well for incubation facilities.

Ken Baker at Scientifix on 1800 007 900 or 03 8540 5900

Email ken@scientifix.com.au

