

Contamination Clean-Up Products



Biocidal ZF™

- Spray Disinfectant for Incubators and Sterile Cabinets in Cell Culture Area.
- Biocidal ZF™ is effective against bacteria, spores and enveloped viruses.
- Surfaces are protected for 2-4 weeks from time of application.
- Active microbiocidal ingredients are non-volatile, preventing contamination of the culture. Cell cultures are protected against contamination and the disinfectant itself.

Biocidal ZF WAK-ZF-1 1 liter **\$230.00 + \$23.00 GST**



AquaClean

- Microbiocidal additive for Water Soluble Heating Bath Fluids.
- Inhibits growth of bacteria, algae, fungi and precipitation of inorganic salts.
- Blue indicator colour. Bacterial contamination is indicated by decolouration of the fluid.
- Does not contain aldehydes or azides.
- Suitable for Stainless Steel – not suitable for baths made of Acrylic, Glass or Aluminium

Aquaclean WAK-AQ-250-50L 250 ml concentrate **\$239.00 + \$23.90 GST**




Acryl AquaClean

- Microbiological Additive for Heating Bath Fluids.
- Compatible with Acrylic, Glass and Aluminium. Universal for all baths.
- Inhibits growth of bacteria, algae and fungi.
- Does not contain aldehydes or azides.

Acryl Aquaclean WAK-AQA-250-50L 250 ml concentrate **\$420.00 + \$42.00 GST**

2012 List Prices*




SystemClean-acidic


SystemClean-acidic and **SystemClean alkaline forte** remove deposits and microbial contamination from analysers and pipettors. Excellent compatibility with materials and fittings used in analytical equipment. **Do not** use for cleaning exchange columns.

Daily Application

SystemClean-alkaline



Weekly Application



CarryOverStop

Microbiocidal Additive for Inter-Assay Wash-Solution. Reduces carry-over effects in analysers and pipettors.

Protolyse enzymatic

Enzymatic Cleaning Agent with Proteolytic Activity for analyzers, pipettors and tubing systems.



Contamination Clean-Up Products

WAK-Chemie Medical GmbH

INNOVATIVE PRODUCTS FOR LIFE SCIENCES



Products for Contamination free Cell Culture include:

Product	Intended target	Effective against	Features
Biocidal ZF™	Spray sterilizer for incubators & sterile cabinets	- Bacteria - Spores - Enveloped viruses	<ul style="list-style-type: none"> • Non-volatile
AquaClean	Sterilizing additive for heating bath fluids	- Bacteria - Algae - Fungi - Precip. of organic salts	<ul style="list-style-type: none"> • Suitable for stainless steel • No aldehydes or azides
Acryl Aqua-Clean	Sterilizing additive for heating bath fluids	- Bacteria - Algae - Fungi	<ul style="list-style-type: none"> • Universal • Suitable for acrylic, glass & aluminium • No aldehydes or azides

Products for Contamination free Analytical Systems include:

Product	Intended target	Removes	Features
SystemClean – Acidic	Analysers & pipettors	- Deposits - Microbial contamination	<ul style="list-style-type: none"> • Aspirate into your system • Excellent for materials and fittings on analytical equipment • Daily application
SystemClean – Alkaline Forte	Analysers & pipettors	- Deposits - Microbial contamination	<ul style="list-style-type: none"> • Aspirate into your system • Effective in inaccessible areas • Daily application
CarryOver Stop	Inter-assay wash solution	- Microbial contamination (Has surfactant properties)	<ul style="list-style-type: none"> • Reduces carry-over effects • Inert to most analytes, enzymes & substrates used in clinical chemistry
Protolyse Enzymatic	Analysers, pipettors, reusable vessels & tubing	- Protein build up	<ul style="list-style-type: none"> • Great for systems in contact with blood, milk or lipoprotein
IronFree	Analysers & equipment	- Iron build up	<ul style="list-style-type: none"> • Compatible with glassware, plasticware & analyser equipment

SystemClean-alkaline forte
SystemClean-acidic
CarryOverStop

WAK-SC-5L
WAK-ASC-5L
WAK-COS-250-5L

5 liters ready to use
250 ml concentrate for 5L
250 ml concentrate for 5L
950 ml stock soln + 50 ml enzyme suspension

inquire
inquire
inquire

Protolyse enzymatic

WAK-PENZ-1L

inquire

IronFree

WAK-IF-1L

1 litre ready to use

inquire

FreeCall: 1800 007 900 info@scientific.com.au

scientific LIFE

12Sfx013

Delivery \$35 + \$3.50 GST. Prices are subject to change without prior notice.