

Mastercycler® X40 PCR Thermocycler

**FREE
REAGENTS**

eppendorf

The New Pulse of Your PCR

Speed, Reliability, and Elegance—Engineered for the Modern Lab

The Eppendorf Mastercycler® X40 is the next evolution in thermal cycling. Designed to be the "daily driver" for busy research and diagnostic labs, it combines the legendary durability of the Mastercycler range with a modern, intuitive interface and sustainable technology. Whether you are performing routine PCR or complex optimization, the X40 delivers results you can trust.

The New Pulse of Your PCR

Exceptional Versatility: The 12-column gradient technology allows you to optimize your PCR protocols in a single run, providing a temperature spread of up to 30°C.

Intuitive VisioNize® Interface

A high-resolution 7-inch touchscreen makes programming effortless. With a "set and forget" workflow, your team can focus on the science, not the hardware.

Superior Block Uniformity

The high-quality aluminum block ensures that every sample across the 96-well plate experiences the exact same thermal conditions, eliminating "edge effects."

Space-Saving Design

With its compact footprint and front-to-back ventilation, the X40 can be placed directly alongside other instruments, maximizing your valuable benchspace.

Sustainability Minded

Features an optimized power-saving mode and eco-friendly packaging, reducing your lab's carbon footprint without sacrificing performance.



MASTER CYCLER X40 – KEY SPECIFICATIONS

| Feature | Specification |
|---------------------|---|
| Block Compatibility | 96-well (0.1 mL or 0.2 mL tubes/plates) |
| Temperature Range | 4°C to 99°C |
| Heating Rate | 3.3°C/s |
| Cooling Rate | 1.5°C/s |
| Gradient Capability | 12-column gradient (up to 30°C) |
| Lid Temperature | 37°C to 110°C |
| Connectivity | Ethernet, USB, VisioNize® Lab Suite |

X40 Mastercycler

#6381000069

List Price
\$6,999 + GST

Limited Time Offer!

Promo Price
\$6,250 + GST



Get a FREE Takara PCR Enzyme + FREE Macherey Nagel Gel & PCR Clean Up Kit

SYNERGY

**Free Samples of
Takara PCR
Polymerase +
Macherey
Nagel Gel & PCR Clean
Up Kit**



Click [here](#) or Scan QR for product details, quotations and demos.



**Free sample of Takara DNA Polymerase (selection to be provided by Scientifix) and Macherey Nagel Gel & PCR Clean Up Kit with the purchase of each X40. Includes Free In Store Delivery. Offer ends June 30th 2026.

PROUD DISTRIBUTOR OF

Clontech Takara cellartis MN MACHEREY-NAGEL



eppendorf

mirus bio

scientifix

1800 007 900 | info@scientifix.com.au | scientifix.com.au