

# Lonza FlashGel™ System



Lonza

## DNA Electrophoresis in Minutes

Get Results in 5 Minutes. No Gels to Prep. No Buffers to Mix.

The Lonza FlashGel™ System is the fastest, most efficient way to separate and recover DNA. By integrating the electrophoresis chamber, transilluminator, and pre-cast gel into one compact unit, it eliminates the traditional bottlenecks of agarose gel electrophoresis.



## Why Choose the FlashGel™ System?

- **Unmatched Speed:** 5 minute separation and recovery. See bands in as little as 2 minutes. Photograph gels at the bench with the optional FlashGel camera.
- **Total Convenience:** The system uses self-contained, pre-cast FlashGel™ Cassettes that include the gel, buffer, and with a highly sensitive stain.
- **Superior Sensitivity:** Detect low levels of DNA (as little as 0.1 ng) with a proprietary stain that is 5–10 times more sensitive than Ethidium Bromide.
- **Real-Time Visualization:** Watch your bands migrate in real-time through the transparent cassette top. The integrated blue-light transilluminator allows for safe viewing without UV damage.
- **Safe Recovery:** The optional FlashGel™ Recovery System allows you to recover DNA samples directly from the cassette, ready for downstream applications—no gel slicing required.

### OFFER 1

**Free Pack of FlashGel DNA Cassettes and FlashGel DNA Marker with FlashGel Device Pack Purchase**

#### FlashGel Device Pack Includes:

- FlashGel Dock: run and visualise samples
- FlashGel Camera: easy image capture
- FlashGel Power Supply: 300V / 400mA
- Windows Software (requires PC not included).

CAT NO.#	Description	Price
57062	FlashGel Device Pack	\$3,922

Select FlashGel DNA Cassette Pack and FlashGel DNA Marker from options on next page.



PROUD DISTRIBUTOR OF



1800 007 900 | [info@scientificx.com.au](mailto:info@scientificx.com.au) | [scientificx.com.au](http://scientificx.com.au)

# Lonza FlashGel™ System



**Lonza**

## OFFER 2

### Free FlashGel DNA Cassettes with FlashGel Power Supply Pack Purchase

#### FlashGel Power Supply Pack Includes:

- FlashGel Dock
- FlashGel Power Supply

CAT NO.#	Description	Price
57069	FlashGel Power Supply Pack	\$2,288



Select FlashGel DNA Cassette Pack from options below.

## OFFER 3

### Free FlashGel Marker with FlashGel DNA Cassettes Purchase

#### Step 1: Purchase Any FlashGel DNA Cassette Below:

Code	Item	No. of Wells	Separation	Pack Size	Price
57023	FlashGel DNA Cassette 1.2%	12+1 well single-tier (13 samples)	50 bp – 4 kb	9 Cassettes	\$359.00
57029	FlashGel DNA Cassette 1.2%	16+1 well double-tier (34 samples)	50 bp – 4 kb	9 Cassettes	\$340.00
57031	FlashGel DNA Cassette 2.2%	12+1 well single-tier (13 samples)	10 bp – 1 kb	9 Cassettes	\$340.00
57032	FlashGel DNA Cassette 2.2%	16+1 well single-tier (34 samples)	10 bp – 1 kb	9 Cassettes	\$340.00



#57023 / 57031




#57029 / 57032

#### Step 2: Select Your Free FlashGel DNA Marker Below:

Code	Item	Size Range	Pack Size	Price
57033	FlashGel DNA Marker	50 bp – 1.5 kb	500 µL	\$418.00
57034	FlashGel DNA Marker	100 bp – 3 kb	500 µL	\$418.00
50473	FlashGel DNA Marker	100 bp – 4 kb	500 µL	\$418.00

Various FlashGel DNA markers run on 1.2% FlashGel Cassette



Click [here](#) or Scan QR for product details, quotations and demos.

**PROMOTION CODE: LFGQ226. Valid until June 30th 2026.**

Please reference promotion code on purchase order. Prices are exclusive of GST. One free FlashGel DNA Cassette Pack and one free FlashGel DNA or RNA marker of choice provided with the purchase of a FlashGel Device Pack. One free FlashGel DNA or RNA marker of choice provided with the purchase of each FlashGel DNA Cassette System. Freight charges apply. For Research Use Only!

PROUD DISTRIBUTOR OF



1800 007 900 | [info@scientificx.com.au](mailto:info@scientificx.com.au) | [scientificx.com.au](http://scientificx.com.au)