

Acryl AquaClean

Microbiological Additive for Heating Bath Fluids
Compatible with Acrylic Glass and Aluminium



Water Baths



Acryl AquaClean inhibits growth of bacteria, algae and fungi.

Acryl AquaClean is long-term compatible with equipment containing acrylic glass and aluminium materials.

Acryl AquaClean does not contain azides or aldehydes.

Application:

Add 5ml AquaClean to 1 liter of heating bath fluid.

Renew fluid and additive every 4 weeks. In case the blue color of fluid has faded, renew earlier.

Compatible for use as an additive in water reservoirs inside cell culture incubators. (Its non-volatile active ingredients cannot invade cell-cultures via air.)

Cat.#	WAK-AQA-250-50L 250 ml concentrate sufficient for 50 liters of fluid	Active Ingredients:	Ketones, Surfactants
Precautions:	Avoid eyes and skin contact	Storage:	Room temperature
Disposal:	Biologically degradable	Materials compatibility:	Universal

scientifix

PO Box 662 South Yarra
VIC 3141 Australia

www.scientifix.com.au

Free call: 1800 007 900
T: +61(0)3 8540 5900
F: +61(0)3 9543 7827
E: info@scientifix.com.au

**WAK - Chemie
Medical GmbH**

