



# Cold Storage

## Preventative Maintenance & Certification Services

### Preventive Maintenance for your Laboratory Fridge / Freezer

- Check electrical safety as per AS 3670 by licensed tester
- Check mains supply and buck/ boost system if fitted
- Check low voltage supplies of control systems and control systems operation
- Carry out any control and display related adjustments necessary
- Check and adjust / repair (if possible) door seals, hinges and locks
- Check and clean condenser
- Check condenser fan / fans
- Check compressor mounts
- Check alarm operation
- Check and replace if necessary alarm back-up battery (specialist batteries supplied at additional cost)
- Check and replace if necessary illumination lamps (specialist lamps at additional cost, eg. fluorescent, UV)
- Check for signs of refrigerant leaks
- Check performance of refrigeration system (where applicable)
- Check heat exchanger temp (where applicable)
- Check compressor temp (where applicable)
- Check condenser super heat (where applicable)
- Check ambient temp and installation

### Certification Services for your Laboratory Fridge / Freezer

- Centre of chamber certification and adjust to specification if necessary / applicable
- Certifications carried out with NATA certified calibrated test equipment and to manufacturer's specifications / procedures / processes

### Quality Service & Support

- All service work carried out by qualified, trained and experienced service staff
- Electrical safety inspections integral to all service work carried out
- Electrical safety certificates / tags available on request for all service work

### Additional Services Available

- Temperature mapping services (please ask for a quote)
- Qualification services available on some makes / models

### **Contact Scientifix and schedule a maintenance & certification service for your cold storage equipment**

Please contact Scientifix for more information, pricing or to request a quote  
For customer service, call 1800 007 900 or email [info@scientifix.com.au](mailto:info@scientifix.com.au)  
More information visit [www.scientifix.com.au](http://www.scientifix.com.au)

**scientifix**