



# Incubators

## Preventative Maintenance & Certification Services

### Preventive Maintenance for your Microbiological Incubator

- Check electrical safety as per AS 3670 by licensed tester
- Check mains supply and max current draw (where appropriate)
- Check condition of electronics and control assemblies, wiring etc
- Check low voltage supplies of control systems and control systems operation
- Carry out any control and display related adjustments necessary
- Check cabinet & chamber fans
- Check door seal, hinges and latches
- Check door switch
- Check water level sensor (if appropriate)
- Check humidification system (if appropriate)
- Check Pressure Equalization Port
- Check chamber access port
- Check all internal and external supply HEPA filters and replace as required (additional cost for filters)
- Check UV lamps and replace if required (parts additional)
- Check for level installation (chamber shelves)
- Check CO<sub>2</sub> detector zero at ambient atmosphere
- Check and record error log and clear – resolve any outstanding issues from error log (as appropriate)
- Check ambient conditions
- Check installation for compliance to manufactures requirements (where appropriate)

### Certification services for your Microbiological Incubator

- Steady state chamber temperature and CO<sub>2</sub> or O<sub>2</sub> measurements to documented procedures (certainty and confidence levels available on request for an additional fee)
- Certifications carried out with NATA certified calibrated test equipment and to manufacturer's specifications / procedures / processes (certainty & confidence level data available on request)

### Quality Service & Support

- All service work carried out by qualified, trained and experienced service staff
- Electrical safety inspections integral to all service work carried out
- Electrical safety certificates available on request for all service work

### Additional Services Available

- Full certification test reports available for all instruments tested, on request, for an additional fee (incl. confidence level and certainty calculations)
- Incubator Qualification services available on request (please ask for a quote)
- Full repair service, depending on parts availability

### Contact Scientifix and schedule a maintenance & certification service for your incubators

Please contact Scientifix for more information, pricing or to request a quote  
 For customer service, call 1800 007 900 or email [info@scientifix.com.au](mailto:info@scientifix.com.au)  
 More information visit [www.scientifix.com.au](http://www.scientifix.com.au)

**scientifix**