



# Total Solution

## Preventative Maintenance & Certification Services

### Comprehensive Maintenance for your Laboratory Instrumentation

- Check electrical safety as per AS 3670 by licensed tester
- Check mains supply and max current draw (where appropriate)
- Check condition of all component assemblies
- Check installation and all service supplies
- Check condition of electronics and control assemblies
- Check low voltage supplies of control systems and control systems operation
- Carry out any control and display related adjustments necessary
- Check and clean refrigeration components (where applicable)
- Check all safety systems
- Check for signs of contamination (where applicable)
- Check performance of all relevant control systems
- Check all cabinet & condenser fan motors
- Check door seal, hinges, stays and latch / interlock
- Check ambient conditions
- Check installation for compliance to manufactures requirements (where appropriate)
- **All repairs, including all parts and labour, for 12months, plus re-certification where appropriate\*, excluding travel and accommodation costs (when required)**  
*\*Please note exclusions may apply as noted on individual contract documents.*

### Certification services for your Laboratory Centrifuges

- Rotational velocity, elapsed time, sample temperature measurements
- Certifications carried out with NATA certified calibrated test equipment and to manufacturer's specifications / procedures / processes\*\*  
 (certainty & confidence level data available on request)  
*\*\* Where appropriate and known*

### Certification services for your Microbiological Incubators

- Steady state chamber temperature and CO<sub>2</sub> or O<sub>2</sub> measurements to documented procedures (certainty and confidence levels available on request for an additional fee)
- Certifications carried out with NATA certified calibrated test equipment and to manufacturer's specifications / procedures / processes (certainty & confidence level data available on request)

### Quality Service & Support

- All service work carried out by qualified, trained and experienced service staff
- Electrical safety inspections integral to all service work carried out, test tags available on request

### Additional Services Available

- Full certification test reports available for all instruments tested, on request, for an additional fee (incl confidence level and certainty calculations)
- GMP/GLP consulting, temp mapping, incubator decontamination services, etc

## Contact Scientifix and schedule maintenance & certification services for your laboratory equipment

Please contact Scientifix for more information, pricing or to request a quote  
 For customer service, call 1800 007 900 or email [info@scientifix.com.au](mailto:info@scientifix.com.au)  
 More information visit [www.scientifix.com.au](http://www.scientifix.com.au)

**scientifix**